


Your Home Therapy Kit.

Acupuncture without Needles



Acupuncture without Needles is a safe, gentle, fast and drug free way to stimulate your body to get out of pain. Adaptive Neuro Stimulation is the key to keep those electronic pulses to continually relax muscles that promotes better blood flow and rapidly releases those all important 'Feel Good' hormones that help you restart the recovery process. Achieving amazing results every time.

Used by professionals, therapists and home users keen on sports or for rehabilitation from chronic health conditions.

At the moment my lower back hurts the most, difficult to reach but with the stretchy strap and circular pad plugged in I was good to go for a twenty minute treatment.

Your home therapy kit comes with a complete set of attachments and a choice of three therapy devices from Denas. Entry level, Intermediate and Advanced devices offers the best flexibility of treatment options making self treatment possible in the comfort of your own home.

I have just received my entry level **Denas 3** device with the full set of attachments and I must say I am impressed with the quality of their new device and equipment.



No stranger to using Adaptive Neuro Stimulation based on the Scenar, RITM technologies, my first device was the AcuScen and I bought my first Denas a couple of years later, back in 2012, all still working today, producing amazing results.

Back then I did not think you could improve on the Denas offerings but they have, its safe, lightweight and a very comfortable fit in the hand.

The Denas 3 is the entry level device for using attachments. All I wanted was to be able to self treat by positioning pads on my back, the area that still hurts the most, I could think about upgrading the device later.

With three modes for treatment, its use is simplicity itself, the frequency of operation is on display along with a timer with the well known healing frequency of 77Hz selected when you switch the device on.

With a choice of three elasticated straps I chose the longest to wrap around my waist and with three body pads to chose from I picked the circular one. With the pad in place and connected to the Denas 3. I turned the device on, selected a twenty minute treatment time and adjusted the intensity to a comfortable level.

In fact I am treating myself as I write this article, its a bit of a luxury to be hands free and receiving a treatment. Moving around is also possible as the circular pad is held comfortably in place with the black elasticated strap.



With this type of technology, with the device on and in contact with your skin you are receiving a whole body treatment with the effects felt immediately and lasting for several days.

The audible feedback is reassuring when you press a button. After turning the Denas 3 on, select the time for a treatment, in this case twenty minutes and then adjust the intensity until you feel the gentle stimulation.

Then I looked at the manual, with clear pictures of the areas to treat, which mode to use and length of treatment time, very easy to understand and implement. It was only then I realised the manual was in Russian and that helped explain some of the symbols on the device. A quick search on the net, English manuals downloaded, all ready to go.

I plan to also use the goggles for five minutes and the foot pad for fifteen minutes with an overall treatment time of forty minutes for this session.



The Denas 3 has three modes of operation, 77Hz, 1040Hz and Med. Pressing the mode button cycles through the modes. In general mode one is the most used, mode two for deep tissue, ligaments, knees, shoulders, elbows, buttocks and mode three is for tonifying or a wellbeing treatment, gets those 'Feel Good' hormones going.

There is also the familiar massage and hair attachments that come as standard in the pack, very effective when you follow the protocols to stimulate the meridians or Acupuncture trigger points.



With this type of technology you can do more if you know Acupuncture, Su Jok, Meridians, Reflexology or researched Bob Becks protocols or Hulda Clark's work.

These pads with the straps give you a convenient way to self treat large areas in the comfort of your own home, hands free.



The two other devices have more options to select from depending on the type of treatment wanted. With the **PCM-4** with 24 preset programmes, ten modes of operation and low frequency selection, the **T3** with ten modes of operation.

When you select a device to use bear in mind what you want to use them for. Maximum flexibility is with the top of the range advanced Denas PCM 4, Spinal treatments, intermediate Denas T3, entry level Denas 3, and the basic Lados.

The Denas 3 is the entry level device for use with attachments.

All these devices come with a 24 month manufacturers warranty. Denas continually improve their technology with their last update in 2014. As you experience and

experiment with what this type of technology can do, you may want to consider an upgrade after a couple of years to their latest devices.

When you order from me at Q + Health you can be assured of prompt delivery with lifetime support by phone, email or skype. As an environmentally conscious company we only issue invoices on request.

If you are lucky enough to live in or visit the Bournemouth and Poole area drop by and see the gadgets, maybe experience how you too can get out of pain.

These really do restart the self recovery process, '[Acupuncture with out Needles](#)'.

<p>Looking for PAIN RELIEF? Is PAIN stopping you TRAIN? Is PAIN putting your life on HOLD?</p>	<p>Accupuncture Without Needles Fast, Effective, Gentle and Drug Free treatment.</p>
	<p>Denas therapy uses Adaptive Neuro Stimulation to relax muscles and rapidly release those Feel Good Hormones. Restarting the self recovery processes.</p>
<p>Home visits now available. For your appointment Call Kevin on 01202 640 264 Today. And experience how quickly you could get out of pain.</p> 	

Pain Relief

Adaptive Neuro Stimulation



Acupuncture without Needles

